

HiSeq Human DNA Resequencing Data Analysis Protocols

Call Variants

文档作者	叶永鑫
提交日期	2012-03-14
更新说明	
文档备注	未完善,需验证

目的: call variants

输入: sorted mapping 结果 bam 文件

库文件: ref.fastn, dbsnp.vcf,

假设已经得到 mapping 结果 bam 文件

(一) 用 samtools、bcftools、vcfutils call variants:

1. filter bam 文件:

留下 unique mapping (对于 bwa 的结果, grep XT:A:U);

滤去 duplicate (Picard-MarkDuplicates);

(PE) 留下 proper paired;

(PE) 留下 insertion size 合理的;(与上面那个有重复?bwa -a 选项?)

(测通的)合并或截短(依赖于建库 fragment 长度,exome seq 有不少测通的 PE reads);

(RNA) 留下链方向正确的

- 1.5. realign? recalibration?
- 2. samtools、bcftools call variants:(http://samtools.sourceforge.net/mpileup.shtml) samtools mpileup -uf ref.fa aln1.bam aln2.bam | bcftools view -bvcg > var.raw.bcf bcftools view var.raw.bcf | vcfutils.pl varFilter -D500 -d 50 > var.flt.vcf # 如果平均深度为 50X samtools call genotype 不准,仅作辅助作用
- (二) 用 GATK call variants:
- 1. MarkDuplicates, Realign, Recalibration(顺序?分 lane 分 sample?)
- 1.1. 先筛出 unique mapping;
- 1.2. Picard-MarkDuplicates

run_picard.sh MarkDuplicates I=read.add_rg.bam O=read.rmdup.bam M=read.rmdup.matrix REMOVE_DUPLICATES=true CREATE_INDEX=true rm read.rmdup.matrix



1.3. GATK-Realign

run_gatk.sh -T RealignerTargetCreator -R ref.fa -I read.rmdup.bam -o read.rmdup.realign.intervals -known dbsnp.vcf run_gatk.sh -T IndelRealigner -R ref.fa -I read.rmdup.bam -targetIntervals read.rmdup.realign.intervals -o read.rmdup.realign.bam rm read.rmdup.realign.intervals

Picard - fix mate information?

1.4. GATK-Recalibration

2. UnifiedGenotyper

(http://www.broadinstitute.org/gsa/gatkdocs/release/org broadinstitute sting gatk walkers g enotyper UnifiedGenotyper.html)

run_gatk.sh -T UnifiedGenotyper -R ref.fa -I read.prepared.bam -glm BOTH --dbsnp dbsnp.vcf -stand_call_conf 30 -stand_emit_conf 30 -o raw.vcf [-alleles dbsnp.vcf] 或

run_gatk.sh -T UnifiedGenotyper -nt 2 -R ref.fa -I read.prepared.bam -D[--dbsnp] dbsnp.vcf -glm SNP -mbq[--min_base_quality_score] 20 -hets[--heterozygosity] 0.001 -I INFO -A[--annotation] AlleleBalance -A DepthOfCoverage -stand_call_conf 30 -stand_emit_conf 10 -dcov 200 -o raw.SNP.vcf run_gatk.sh -T UnifiedGenotyper -nt 2 -R ref.fa -I read.prepared.bam -D dbsnp.vcf -glm INDEL -mbq 20 -indelHeterozygosity 0.000125 -I INFO -A AlleleBalance -A DepthOfCoverage -stand_call_conf 30 -stand_emit_conf 10 -dcov 200 -o raw.indel.vcf (一起 call 较快)

或

run_gatk.sh -T UnifiedGenotyper -R ref.fa -I read.prepared.bam -glm BOTH
[-B:alleles,VCF dbsnp.vcf -BTI alleles] -B:dbsnp,VCF dbsnp.vcf -stand_call_conf 50
-stand_emit_conf 10 -dcov 1000 --min_base_quality_score 30 -A DepthOfCoverage -A
AlleleBalance -o raw.vcf -metrics raw.metrics

3. Variant Select

run_gatk.sh -T SelectVariants -R ref.fa --variant raw.vcf -selectType SNP -selectType MNP -o raw.snp.vcf run_gatk.sh -T SelectVariants -R ref.fa --variant raw.vcf -selectType INDEL -o raw.indel.vcf



4. Variant Filtration

run_gatk.sh -T VariantFiltration --clusterWindowSize 10 --filterExpression
"MQ0>=4&&((MQ0/(1.0*DP))>0.1)" --filterName "HARD_TO_VALIDATE" --filterExpression
"DP<10" --filterName "LowCoverage" --filterExpression "QUAL<30.0" --filterName
"VeryLowQual" --filterExpression "QUAL>30.0&&QUAL<50.0" --filterName "LowQual"
--filterExpression "QD<5.0" --filterName "LowQD" --filterExpression "SB>-0.10" --filterName
"StrandBias" -B:mask,VCF raw.indel.vcf --maskExtension 0 --maskName Indel -R ref.fa
-B:variant,VCF recal.snp.vcf -o flt.snp.vcf
(対 exome seq 不推荐 Qual 过滤)
或

(三) Report

Ti、Tv

VariantEvalu

QC countLoci countPairs ...

GATK-Phasing,对于家系重要 GATK-somaticIndelDetector dIndel?